DEMSMD Error Correction 2024-01-23:

1. 9120 Epinephrine – Under the adult epinephrine infusion the label concentration is incorrect

Adult Infusion:

- Mix: inject 1 mg epinephrine into 1000 mL Normal Saline bag to achieve 1 mcg/mL concentration (This means 1 mL of 1mg/mL [1:1000] or 10 mL of 0.1mg/mL [1:10,000] either way 1 mg of drug). Use macro drip set.
- Adult IV/IO: Begin IV/IO infusion wide open to gravity to give small aliquots of fluid.
 Typical volumes are less than 100 mL of total fluid, as typical doses are expected to be < 100 mcg. Titrate to desired hemodynamic effect with goal BP of > 90 mmHg systolic, improved respiratory status (bronchodilation), and improved perfusion/mentation.
- Label: Apply label to bag noting epinephrine 10mcg/mL

Protocol updated to read label adult infusion saline bag as 1 mcg/mL

Adult Infusion:

- Mix: inject 1 mg epinephrine into 1000 mL Normal Saline bag to achieve 1 mcg/mL concentration (This means 1 mL of 1mg/mL [1:1000] or 10 mL of 0.1mg/mL [1:10,000] either way 1 mg of drug). Use macro drip set.
- Adult IV/IO: Begin IV/IO infusion wide open to gravity to give small aliquots of fluid.
 Typical volumes are less than 100 mL of total fluid, as typical doses are expected to be < 100 mcg. Titrate to desired hemodynamic effect with goal BP of > 90 mmHg systolic, improved respiratory status (bronchodilation), and improved perfusion/mentation.
- Label: Apply label to bag noting epinephrine 1 mcg/mL